

Dehydrated Skin 2 step fixes

Dehydrated Skin 2 step fixes

[Dehydrated Skin Products](#)

Dry and Dull Dehydrated Skin

Apply 1-2 pumps of [Citrus Stem Cells + Hyaluronic Acid](#)

Apply 1-2 pumps [Cell Restore + Phytic Acid Moisturizer](#)

Use AM & PM

Rough Texture Dehydrated Skin

Apply 1 pump [Hydrating Oil w/Vitamin C](#)

Apply 2 [pumps Antioxidant Rich SPF 30](#)

Use in AM

Brighten and Normalize Dehydrated Skin

Apply 2 pumps [Bakuchoil + Niacinamide + Squalane Serum](#)

Apply 2 pumps Tri-Peptide with Vitamin C

Apply PM

Brighten & Smooth Dull, lackluster, Tired, Dehydrated Skin AM

Apply 3 drops of [Complex Vit. C + Antioxidants Serum](#)

Layer a dime size of [Moisture Lock](#)

Apply at AM

Brighten & Smooth Dull, Lacluster, Tired, Dehydrated Skin PM

Apply 3 drops of [Hyaluronic Acid + Antioxidants](#)

[Layer Detoxifying 3in1 Hydrating Masque](#)

Apply at PM